
Problem: Calibration (easyTouch Display)

- Calibration of the touchscreen is lost (e.g. after changing the date/time, tones, transfer)
- Calibration is not applied after recalibration/restart

Cause

- The device stops working while storing to the registry
- The device is switched off while storing to the registry



The next time you start, you will be prompted to perform a recalibration.

Action

➤ Calibration

1. A cross appears at the center of the touchscreen.
2. Press the center of the cross (in total 5 times).
3. Follow the rest of the crosses and press the center of the cross each time. When doing so, do not touch any other place on the touchscreen.
4. Confirm the action by pressing the touchscreen anywhere within 30 seconds after completion of the calibration.

➤ Saving the settings to the registry

1. Open the 'Settings' page.
2. Open the 'Transfer' menu.
3. Press the 'Network' button.
4. Save the IP address.
5. Save the settings.



**Wait until the saving process is finished.
Never switch off, before saving process has finished completely.**